

124mm

265mm

Dr. Odin®
DIGITAL THERMOMETER
Instruction Manual

Model No.: ODT311

Intended Use: It is a device used to record body temperature.

Warning :

- ⚠ Read instructions thoroughly before using digital thermometer.
- ⚠ Choking Hazard: Thermometer cap and battery may be fatal if swallowed. Do not allow children to use this device without parental supervision.
- ⚠ Do not use thermometer in ear. Designed use is for oral, rectal, and armpit (axilla) readings only.
- ⚠ Do not place thermometer battery near extreme heat as it may explode.
- ⚠ Note: Use of the probe cover may result in a 0.1°C (0.2°F) discrepancy from actual temperature.
- ⚠ It is recommended the performance should be checked every two years.
- ⚠ Remove battery from the device when not in operation for a long time.
- ⚠ The use of temperature readings for self-diagnosis is dangerous. Consult your doctor for the interpretation of results. Self-diagnosis may lead to the worsening of existing disease conditions.
- ⚠ Do not attempt measurements when the thermometer is wet as inaccurate readings may result.
- ⚠ Do not bite the thermometer. Doing so may lead to breakage and/or injury.
- ⚠ Do not attempt to disassemble or repair the thermometer. Doing so may result in inaccurate readings.
- ⚠ After each use, disinfect the thermometer especially in case the device is used by more than one person.
- ⚠ Do not force the thermometer into the rectum. Stop insertion and abort the measurement when pain is present. Failure to do so may lead to injury.
- ⚠ Do not use thermometer orally after being used rectally.
- ⚠ For children who are two years old or younger, please do not use the devices orally.

PLEASE READ CAREFULLY BEFORE USING


This digital thermometer provides a quick and highly accurate reading of an individual's body temperature. To better understand its functions and to provide years of dependable results, please read all instructions first.

This thermometer is manufactured according to the system of EN12470-3:2000+A1:2009 ASTM E1112-00(Reapproval 2011). And the manufacturer is EN ISO 9001:2008 and EN ISO 13485:2012/AC:2012.

CONTENTS

1 Thermometer, 1 Owners Manual, 1 Storage Case

PRODUCT ILLUSTRATION



PRECAUTION



- * The performance of the device may be degraded should one or more of the following occur:
 - Operation outside the manufacturer's stated temperature and humidity range.
 - Storage outside the manufacturer's stated temperature and humidity range.
 - Mechanical shock (for example, drop test).
 - Patient temperature is below ambient temperature.
- * Portable and mobile RF communications such as microwave ovens and cell phones, which generate strong electrical or electromagnetic fields, near the device, as they may become interference source to the device.

SYMBOL EXPLANATION

---	Direct Current	LOT	Batch code
⊠	Type BF Applied Part	⏏	Manufacturer
⚠	Caution, consult accompanying documents		

SPECIFICATIONS

Type:	Digital Thermometer (Not Predictive)
Range:	32.0°C, 42.9°C (90.0°F - 109.9°F) (°C/°F chosen by manufacturer)
Accuracy:	±0.1°C, 35.5°C, 42.0°C (±0.2°F, 95.9°F - 107.6°F) ±0.2°C under 35.5°C or over 42.0°C (±0.4°F under 95.9°F or over 107.6°F) at ambient temperature of 18°C to 28°C (64.4°F to 82.4°F)
Display:	Liquid crystal display, 3 digits
Memory:	For storing the last measured value
Battery:	One 1.5 V DC button battery (size LR41 or SR41, UCC 392)
Battery life:	Approx 200 hours
Dimension:	12.4cm x 1.8cm x 0.9cm (L x W x H)
Weight:	Approx. 10 grams including battery
Using environment:	Temperature: 5°C - 40°C (41°F - 104°F) Relative humidity: 15%-95%RH
Storage and transportation condition:	Temperature: -20°C ~ 55°C (-4°F - 131°F) Relative humidity: 15%-95%RH
Ingress Protection Rating:	IP 22
Classification:	Type BF ⊠

DIRECTIONS

1. Press the *On/Off* Button next to LCD display. A tone will sound as the screen shows **1888.5** followed by last recorded temperature. After showing the self-test temperature, the thermometer is now in the testing mode.

*Note: If the measured temperature is less than 32°C or 90°F, the LCD will display Lo, if the measured temperature is higher than 42.9°C or 109.9°F, the LCD will display Hi.

2. Position thermometer in desired location (mouth, rectum, or armpit.)
 - a) **Oral Use:** Place thermometer under tongue as indicated by "✓" position shown in Figure 2. Close your mouth and breathe evenly through the nose to prevent the measurement from being influenced by inhaled/exhaled air. Normal temperature between 35.7°C and 37.3°C (96.3°F and 99.1°F)
 - b) **Rectal Use:** Lubricate silver probe tip with petroleum jelly for easy insertion. Gently insert sensor approximately 1cm (less than 1/2") into rectum. Normal temperature between 36.2°C and 37.7°C (97.2°F and 99.9°F).
 - c) **Armpit Use:** Wipe armpit dry. Place probe in armpit and keep arm pressed firmly at side. From a medical viewpoint, this method will always provide inaccurate readings, and should not be used if precise measurements are required. Normal temperature between 35.2°C and 36.7°C (95.4°F and 98.1°F).




Figure 2

3. The degree sign flashes throughout the testing process. When flashing stops an alarm will beep for approximately 10 seconds. The measured reading will appear on the LCD simultaneously. The minimum measurement time until the signaling tone (beep) must be maintained without exception. The measurement continues even after the buzzer notification. The displayed temperature will not change when thermometer is removed from its testing position.

*Note: Normally the buzzes are "Bi-Bi-Bi-Bi-"; Alarm beeps more rapidly when temperature reaches 37.8°C (100°F) or higher, and the buzzes are "Bi-Bi-Bi---- Bi-Bi-Bi----- Bi-Bi-Bi"

4. To prolong battery life, press the *On/Off* Button to turn unit off after testing is complete. If no action is taken, the unit will automatically shut off after around 10 minutes.
5. Store the thermometer in its protective case.

BATTERY REPLACEMENT

1. Replace battery when "⏏" appears in the lower right corner of LCD display.
2. Pull battery cover off as shown in Figure 3.
3. Gently pull out plastic circuit board with battery chamber approximately 1 cm (slightly less than 1/2") (See Figure 4)
4. Use pointed object such as a pen to remove old battery. Discard battery lawfully. Replace with new 1.5V DC button type LR41 or SR41, UCC392, or equivalent. Be sure battery is installed with "+" polarity facing up. (See Figure 5)
5. Slide battery chamber back into place and attach cover.




Figure 3




Figure 4




Figure 5

CLEANING AND DISINFECTION

- Wipe the thermometer with a soft clean cloth.
- For stubborn stains, wipe the thermometer with a cloth that has been dampened with water or a neutral detergent solution and then wring thoroughly. Finish by wiping with a soft dry cloth.
- Observe the following to prevent damage to the thermometer.
 - Do not use benzene, thinner, gasoline or other strong solvents to clean the thermometer.
 - Do not soak the sensing section in alcohol for long periods of time or attempt to sterilize it using hot water (water at a temperature of 50 °C (122 °F) or higher).
 - Do not use ultrasonic washing to clean the thermometer.

Passim Lifesciences Ltd. warrants this units free of any defect in material, workmanship and operation, under normal use, for a period of one year from the date of purchase. If any part become defective during the warranty period Passim Lifesciences Ltd. will repair or replace (if not repairable) the same, free of cost. This shall not apply to any parts that are considered as expandible or deteriorable in the course of normal use. Passim Lifesciences Ltd. shall be relieved of any liability and warranty shall cease to apply if :

- This is not used in accordance with the instructions in the operational manual.
- It is used with any equipment not complying with the specification of this unit.
- It is not regularly maintained.
- The unit is disassembled, repaired or operated by person not authorized by Passim Lifesciences Ltd.
- Damage caused due to negligency.
- The unit is operated in corrosive materials or in the harmful atmosphere.
- The warranty card is not filled completely and produced at the time of warranty claim.

Marketed By:
Odin Healthcare Pvt. Ltd.
Plot No. 45, Ind. Area, Phase-II,
Panchkula - 134113 (HR) - INDIA
For any Complaint/ Suggestion
please contact :
Customer Care No.: 1800 309 3009
(Timing: 9am - 7pm, Mon. - Sat.)
Email ID: customercare@drodin.in
Website: www.drodin.in

WARRANTY CARD	
Model No.	Lot No.
Invoice No.	Date of Purchase.
Purchased By:	Contact No.:
Address:	
Dealer's Name:	
Dealer's Sign & Stamp	

* Terms & Conditions apply